


### **Analyze Particles**

#### **Image Processing & Analysis for Life Scientists**

Olivier Burri, Romain Guiet & Arne Seitz









### **Analyze Particles**





- Using Analyze Particles
- Displaying Results
- Output to Roi Manager

Hi! In this video, we will cover the analyze particles command, and some of its possible outputs. We will focus on how to use the particle analyzer in different ways. How to display results. And where they come from. As well as introduce the ROI manager.


- Summary -

	回海銀河
5	22000
=	
5	<b>间数44</b> 4

### **Analyze Particles**








Now, you have already seen these images before. These steps show us how to get from an image to objects. And we want to cover that last part down there. But to begin, let's apply some of the knowledge we've learned over the previous weeks to improve the quality of our detection. We've seen how filtering can help us remove noise. So in this case, let's apply a median filter to the initial image. Morphological operations in binary images can help us separate the particles further, like the watershed filter that we apply here. And this will lead to our object detection being a lot cleaner. So starting with this mask, let's see how we can get objects using the analyze particle command.

	- Summary ————————————————————————————————————
	Julilliary
国 成	

5.4. Particle Analysis 3 of 16

# In Fiji





The command can be found under the Analyze menu. And the window looks like this. This command is closely linked to the set measurements command as we will see later on. There are many options here.

- Notes

Summary —			
		- Summary	
		Carrinary	
	国際教練国		
\$297.6236 6.847.6346	(本)		
NS#接種			
	ASANSESE		
	国外投资数		

## In Fiji





For the purpose of this video, we will only focus on the size and the circularity fields that allow us to filter out particles that would or would not meet this criteria. We will cover a few outputs of the Show field, and touch on the Display Results, and Add to Manager check boxes.


Notes -	

	— Summary —————	
The Marie To		

5.4. Particle Analysis 5 of 16








The size filter lets a select particle to keep within a certain range. Notice that because our test image is calibrated, the size is shown in calibrated units.



	Summary —
	Summary






So if we only want to keep bacteria above 1 micron, we can do so. To see the results we set the Show field to Masks. This will create a neaw mask image with only the particles that match the size criteria, ready for further processing.

- Summary —	

5.4. Particle Analysis 7 of 16





Similarly, we can chose to only keep bacteria that are round by setting the allowed circularity to be between 0.8 an 1. To get more information on this measurement types, you can see the documentation table on the ImageJ the website. And the link is available below the video.


Summary —	





Input





Output **Glasbey Inverted** LUT



To get the results we showed at the beginning, the option is called Count Masks under the 'Show' field. This will give you an image with each individual particle, numbered with a unique number. And then you can color it using the appropriate color table. In this case we use the inverted Glasbey table.


Notes —

	- Summary	
	Summary	
ELENATE		
<b>国場織国</b>		
<b>光彩的数形</b> 型		

9 of 16 5.4. Particle Analysis

#### **Analyze Particles: Overlay**






So far, we did not see any data coming out of this analysis. So, to do that we can use the Display results check box to obtain information on each detected particle. A very useful feature to check the results is to use the show overlay mode. We run Analyze particles and then we select Show Overlay. When we click 'Ok' we get the results table. And if we use the rectangle tool whatever particle we select on our image will give us the corresponding line in the results table. Unfortunately, it doesn't work the other way. You cannot click on a result and then see what particle it corresponds to.

N	lotes ————
1	.0.00

- Summary

## **Analyze Particles: Overlay**





If you wish to remove the overlay simply go to image overlay remove overlay.


10103

1		Summary
	■統治(国 320年66年	
1		

5.4. Particle Analysis 11 of 16

### **Analyze Particles: Results**





The current results table has 2 columns. So let's see where these are coming from. Like we said, the set measurements command dictates what gets measured when analyze particles is run.




Notes

Guilliary

### **Analyze Particles: Results**





We see that right now only the Area and Display label are selected. We recommand you always have Display label cheked to make sure you know what image is being analyzed. How to select and manipulate the results will be covered in a subsequent video. But for now, it's good that you note how to obtain it with the analyze particle command.


[ NOIES		
_		
-		
_		
_		
$oxed{}$		

	Summary —	
<b>国政的经验</b>		

5.4. Particle Analysis 13 of 16

### **Analyze Particles: ROI Manager**






The last option is with regards to sending the detected particles to something called the ROI Manager. This will create ROI (regions of interest) which are overlaid on the original images. These can then be used to visually inspect each particle, and see for example if the detection was correct. Furthermore, these regions are now independent of the masks image and can be to re-applied to other images, or other channels of the same image for example.

1	
1	
1	
1	
I	
1	
I	
I	
1	
1	
I	
1	
I	
1	
l	
l	
1	
1	
1	
1	
I	
1	
I	
1	
1	
$\overline{}$	

- Summary —	

# **Analyze Particles: ROI Manager**





ROI will be covered more extensively later. So don't worry if it seems that we've skipped through it.

110100

Summary —	

5.4. Particle Analysis 15 of 16

#### **Conclusion**





- Using Analyze Particles
- Displaying Results
- Output to ROI Manager

Alright, that's the end for this video. We saw how to analyze the results for a binary mask. And how to perform filtering based on size and circularity. As well as what kind of outputs are available from the analyze particles command. We touch base on the concept of Results tables and Region of Interest. Both of which can be produced by this command. Goodbye! And enjoy the exercises!

Summary

Notes

	同語なる同
)	Lette -Return 1
Н	POUR TRANSPORT
П	A COMPANY OF THE
.	
: 1	
:	
	E