


### **Results Tables**

#### **Image Processing & Analysis for Life Scientists**

Olivier Burri, Romain Guiet & Arne Seitz









### Outlook





- Setting Measurements
- Getting/Setting Results
- Saving Results

Hi! And welcome to this video on results tables. We will cover how to select a measurement we want to display into our results tables. We will also see how to get results from a table. And how to build our own results columns using the macro language. Finally, we will see in which format results tables are stored when they are saved.

-	

- Summary

6. ROIs & Results 2 of 13

# From Images to Data





Image Filtered

Mask Watershed

Object

Over the previous weeks, we have reached the point where we now have objects that represent what we're interested in measuring.

Notes

	———— Summary ————	
1255		


6.3. Results Tables 3 of 13

## From Images to Data





Object



bacteria\_crop-3.tif 0.604 bacteria\_crop-3.tif 2.000 bacteria\_crop-3.tif

hacteria cron-3 tif

1.396

1.856

Measurements

So, for example, in the case of bacteria here, say, we would be interested in measuring their size! So now we need to go to the extra steps and see how and where we can extract this information from these objects.

 Notes —
Notes —


	— Summary ———	
	Carrinary	
回鉄返回		

6. ROIs & Results 4 of 13

# **Setting Measurements**








We saw the analyze particle tools, and the ROI manager in the previous weeks. And we saw that these can output results tables. We also mentioned how these are linked to the set measurements command, like the label and area here.

	Notes —
'	10100

6.3. Results Tables 5 of 13

## **Setting Measurements**





Activating the check boxes simply create new columns in the results table that will host all the measurements for us.

	Notes
l	


- Summary	

6. ROIs & Results 6 of 13

## **Setting Measurements**







Note that if you already ran measurements, checking new boxes will not recompute the measurements. You're going to have to run either analyze particles or the ROI Manager measure command again in order to populate those columns.



Notes

	Summary —	
	Sullillary	
0.00 (0.00 ft)   0.00 (0.00 ft)   0.00 (0.00 ft)		


6.3. Results Tables 7 of 13


### **Managing Results**





- If Results Table Exists
  - Appends to Results
  - Close or Clear Results to get new table





Now, as a default, results in Fiji will always be appended at the end of an existing table. That is, new results will be creating new rows into the results table that is already opened. However you can use the clear results command from the analyze particles menu, or close the results table, or use clear results, if you need a fresh new empty table. Now what if we need more results or more data? Or what if we want to compute a new result based on some existing column?

Notes —

- Summary —	

6. ROIs & Results 8 of 13

#### What if we Need More?






#### Macro Language Results Functions

Example: Binning by size

value = getResult("Column Name", row);
setResult("Column Name", row, value);



For example, what if we want to classify bacteria into 3 categories: Small, medium and large based on how large their areas are? In this case, the macro language has a series of functions available that allows to navigate, modify, and create new results. These are all available in the image macro functions page. But the most important ones are 2 functions called Get Results and Set Results. Get results will turn the value of the column defined by the name, and at the selected row. Set results normally will add a value to a given column, at a given row, or create a new column if it doesn't exist yet. So if we take the results of the area column, and we choose a label based on this value (say small it it's smaller than 1 micron medium, if it's between 1 and 2 microns and large if it's above 2) we can create a new measurement and add it to the results table. You can see the marco code for this example at the end of the video.



The 258

6.3. Results Tables 9 of 13

Summarv

# **Getting Statistics**





- Fiji is not a Data Analysis Tool
  - Use your favorite software











However, we will not recommend you use Fiji for analyzing your data, or getting statistics. The sofware is useful for image analysis and for the extraction of results from images. Any kind of further analysis, specially statistical analysis, should be carried out in a specialized sofware, like Excell R Igor, just to name a few.

Notes

	Summary —
)	
1	

6. ROIs & Results 10 of 13


# **Saving Results**





#### Saved as CSV Files

• Can be easily imported into any software



But for that, we need to be able to export the results from Fiji. So, this is very easily done via the file Saved as command. This will create a common seperated file wherever you choose. Most software packages can import CSV Files very easily.

Notes

- Summary —	

6.3. Results Tables 11 of 13


## **Saving Results**





#### Saved as CSV Files

- · Can be easily imported into any software
- Table Options: Edit > Options > Input/Output



Now if you are not getting the same results as the examples that I've shown here, check the input output option of Fiji under the Edit option input output menu. Here you can select the extension for which you want the results tables to be saved. And whether the column names and the row numbers should be saved along side with your data.

L	
ľ	
ŀ	
ľ	
ŀ	
ľ	
ŀ	
ľ	
ŀ	
ľ	
ŀ	
ľ	
ŀ	
ľ	
ŀ	
ľ	
(	

Notes

	Summary	
ļ		

6. ROIs & Results 12 of 13

### **Conclusion**





- How to Set Measurements
  - · Need to run measurement after setting
- Getting/Setting Results
  - getResult() setResult()
- Saving Results
  - Result Options

In this video, we've covered how to get measurements out of an image. We've seen that there is a little detail in which you have to be careful to re-run the measurements, if you select a new one from the set measurements command. We've also seen how to get and set results from column tables, which allows us to create our own custom results. And we've seen that results can be saved in a convenient CSV format. See you next time!

Notes

	- Summary	
	Sullillary	
Total (		

6.3. Results Tables 13 of 13